



**Bio Optics**  
Product Catalog

The logo features a stylized blue swirl with five colorful petals (yellow, pink, purple, teal, and light blue) radiating from the top. To the right of the swirl, the text "Bio Optics" is written in a large, bold, dark blue sans-serif font, with "Product Catalog" in a smaller, orange sans-serif font below it.

*"New runners of Medical market"*

# BIO Plug™

(Dissolvable Lacrimal Punctum Plug)



## BIO Plug

- Use to Polydioxanone Material
- Violet Color
- Never touch the eye
- **Comfortable** after proper placement
- Absorption of plug approximately 180days
- Available in two size:0.4mm and 0.5mm diameter

### [ Indication ]

- Dry eye
- Refractive surgery
- Contact lens wear
- Seasonal allergy
- Cataract Surgery



0.4mm or 0.5mm



10 pcs/box



# BIO Intubation Set™

(Lacrimal Intubation Set)


## BIO Lacrimal Intubation Set

- New state art of Lacrimal Intubation system.
- Bio compatible silicone tubing.

5 sets/box


### Stainless steel probes with silicone tubing.

BS-1001



Probes – 11 cm length X 0.40 mm diameter (4 1/4 in X 27G)  
 Silicone tubing – 30 cm length, 0.30 mm inside diameter,  
 0.64mm outside diameter

BS-1002



Probes – 11 cm length X 0.60 mm diameter (23G x 7 in)  
 Silicone tubing – 30 cm length, 0.30 mm inside diameter,  
 0.64 mm outside diameter

## BIO Lacrimal Intubation Set

- New Design of Lacrimal Intubation probe
- Easy to Use.

5 sets/box

### Stainless steel probes with PDO tubing.

BP-2001



Probes – 11 cm length X 0.40 mm diameter (4 1/4 in X 27G)  
 PDO tubing – 35 cm length, 0.30 mm inside diameter,  
 0.64mm outside diameter

BP-2002



Probes – 11 cm length X 0.60 mm diameter (23G x 7 in)  
 PDO tubing – 35 cm length, 0.30 mm inside diameter,  
 0.64 mm outside diameter

# BIO Color Tear Test™

(Schirmer Strip)

- Easy to use
- Blue-black dye line
- Different shapes for right(R) and left(L) eye



100 pcs/box




## Critical Value

- Normal which is 10mm–30mm wetting of the paper after 5minute.
- Abnormal which is < 10mm wetting of the paper after 5minute.

# BIO PVA Spear™

- PVA Eye Spear
- Stiff without water and soft with water
- Provide high safety to the cornea
- Micro point absorbent
- Lint free
- Plastic handle



100 pcs/box

# BIO Micro Syringe Filter™

- Bio PES Syringe Filter
- Specification
  - 0.2µm x 25mm
  - Safety to PH1~14
  - Hydrate Membrane Filter
  - Gamma Radiation



50 pcs/box



**Bio Optics Co., Ltd.**

No.1003, SICOX Tower,484 Dunchon-daero, Jungwon-gu, Seongnam-si,  
Gyeonggi-do, KOREA 13229  
T.82-31-747-4850 F.82-31-747-4820

**"New runners of Medical market"**